

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\Localizer
 TA:0:14 PAT:Off Voxel size:0.5x0.5x7.0 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	On
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	L0.0 A20.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	7.0 mm
TR	8.6 ms
TE	4.00 ms
Averages	2
Concatenations	3
Filter	Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	91 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
TD	0 ms
Prescan Normalize	Off
Normalize	On
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	L0.0 A20.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SRFExcit 1H	43.920 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	2D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	320 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Liver registration	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\AAHead_Scout_32ch-head-coil
TA:0:14 PAT:3 Voxel size:1.6×1.6×1.6 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP
AutoAlign	Head

Contrast

Flip angle	8 deg
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude

Resolution

Base resolution	160
Phase resolution	100 %
Phase partial Fourier	6/8
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	69 %
Slice partial Fourier	6/8

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
Multi-slice mode	Sequential
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SRFExcit 1H	41.037 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Weak
Contrasts	1
Bandwidth	540 Hz/Px
RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Time to center	6.2 s
Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Contrasts	1
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\Localizer_aligned
 TA:1:00 PAT:Off Voxel size:0.5×0.5×10.0 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	On
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	3
Slices	9
Dist. factor	80 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	8.6 ms
TE	4.00 ms
Averages	2
Concatenations	15
Filter	Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	91 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
TD	0 ms
Prescan Normalize	Off
Normalize	On
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Nr. of slice groups	3
Slices	9
Dist. factor	80 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SRFExcit 1H	43.920 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Dimension	2D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	320 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Liver registration	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\BIAS_BC
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	280.0 ms
TE	1.02 ms
Averages	1
Concatenations	1
Filter	Normalize
Coil elements	BC

Contrast

Magn. preparation	None
Flip angle	3 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	On
Cut off	20
Width	4
Unfiltered images	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SRFExcit 1H	46.167 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3.6 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	BC
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\BIAS_32CH
 TA:0:28 PAT:Off Voxel size:2.0×2.0×2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	280.0 ms
TE	1.02 ms
Averages	1
Concatenations	1
Filter	Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	None
Flip angle	3 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	On
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SRFExcit 1H	46.167 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3.6 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\SpinEchoFieldMap_PA
 TA:0:31 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	7700 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
AddCSaCSatNS 1H	68.283 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\SpinEchoFieldMap_AP
 TA:0:31 PAT:Off Voxel size:2.0×2.0×2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	7700 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
AddCSaCSatNS 1H	68.283 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\rfMRI_REST_AP
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	720 ms
TE	37.00 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
MBExc 1H	504.090 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	20
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\rfMRI_REST_PA
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	720 ms
TE	37.00 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
MBExc 1H	504.090 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	20
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\T1w_MPR
 TA:6:38 PAT:2 Voxel size:0.8×0.8×0.8 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 P3.0 H6.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	23.1 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	2400.0 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1060 ms
Flip angle	8 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	23.1 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	0.00 deg
F >> H	256 mm
A >> P	240 mm
R >> L	167 mm
Frequency 1H	123.173991 MHz
Correction factor	1
ExcitWEns 0 1H	61.556 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1060 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	210 Hz/Px
Flow comp.	No
Echo spacing	8.1 ms
Turbo factor	256
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
MapIt	None
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\T2w_SPC
 TA:5:57 PAT:2 Voxel size:0.8×0.8×0.8 mm Rel. SNR:1.00 :spc

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Position	L0.0 P3.0 H6.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	3200 ms
TE	564.0 ms
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	None
Restore magn.	Off
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Allowed
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	208
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off
Restore magn.	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.173991 MHz
Correction factor	1
SPC_nsExc 1H	230.833 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Trigger delay	0 ms
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Bandwidth	744 Hz/Px
Flow comp.	No
Allowed delay	0 s
Echo spacing	3.86 ms
Adiabatic-mode	Off
Turbo factor	314
Echo train duration	1166 ms
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Organ under exam.	Standard
Tissue T1	940 ms
Tissue T2	100 ms

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\DWI_dir98_AP
TA:5:35 PAT:Off Voxel size:1.5×1.5×1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	3222 ms
TE	89.20 ms
Averages	1
Multi-band accel. factor	4
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	210 mm
A >> P	210 mm
F >> H	138 mm
Frequency 1H	123.173991 MHz
Correction factor	1
ExtExciteRF 1H	119.065 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Free echo spacing	Off
Echo spacing	0.69 ms
SIR accel. factor	1
EPI factor	140
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\DWI_dir98_PA
 TA:5:35 PAT:Off Voxel size:1.5×1.5×1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	3222 ms
TE	89.20 ms
Averages	1
Multi-band accel. factor	4
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	210 mm
A >> P	210 mm
F >> H	138 mm
Frequency 1H	123.173991 MHz
Correction factor	1
ExtExciteRF 1H	119.065 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Free echo spacing	Off
Echo spacing	0.69 ms
SIR accel. factor	1
EPI factor	140
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\SpinEchoFieldMap_PA
 TA:0:31 PAT:Off Voxel size:2.0×2.0×2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	7700 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
AddCSaCSatNS 1H	68.283 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\SpinEchoFieldMap_AP
 TA:0:31 PAT:Off Voxel size:2.0×2.0×2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	7700 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
AddCSaCSatNS 1H	68.283 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\rfMRI_REST_AP
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	720 ms
TE	37.00 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
MBExc 1H	504.090 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	20
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\rfMRI_REST_PA
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	720 ms
TE	37.00 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	208 mm
A >> P	208 mm
F >> H	144 mm
Frequency 1H	123.173991 MHz
Correction factor	1
MBExc 1H	504.090 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.58 ms
SIR accel. factor	1
EPI factor	104
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	20
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\DWI_dir99_AP
 TA:5:38 PAT:Off Voxel size:1.5×1.5×1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	3222 ms
TE	89.20 ms
Averages	1
Multi-band accel. factor	4
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	210 mm
A >> P	210 mm
F >> H	138 mm
Frequency 1H	123.173991 MHz
Correction factor	1
ExtExciteRF 1H	119.065 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Free echo spacing	Off
Echo spacing	0.69 ms
SIR accel. factor	1
EPI factor	140
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\HCP-Lifespan\Lifespan\Session 1\DWI_dir99_PA
 TA:5:38 PAT:Off Voxel size:1.5×1.5×1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	3222 ms
TE	89.20 ms
Averages	1
Multi-band accel. factor	4
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off

Geometry

Nr. of slice groups	1
Slices	92
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	210 mm
A >> P	210 mm
F >> H	138 mm
Frequency 1H	123.173991 MHz
Correction factor	1
ExtExciteRF 1H	119.065 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Free echo spacing	Off
Echo spacing	0.69 ms
SIR accel. factor	1
EPI factor	140
Gradient mode	Performance
RF spoiling	Off
Online multi-band recon.	Online
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	2
b-value 1	0 s/mm ²
Diff. weighted images	On
Trace weighted images	On
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D**Table of contents**

```
\\USER
  | HCP-Lifespan
  |   | Lifespan
  |   |   | Session 1
  |   |   |   | Localizer
  |   |   |   | AAHead_Scout_32ch-head-coil
  |   |   |   | Localizer_aligned
  |   |   |   | BIAS_BC
  |   |   |   | BIAS_32CH
  |   |   |   | SpinEchoFieldMap_PA
  |   |   |   | SpinEchoFieldMap_AP
  |   |   |   | rfMRI_REST_AP
  |   |   |   | rfMRI_REST_PA
  |   |   |   | T1w_MPR
  |   |   |   | T2w_SPC
  |   |   |   | DWI_dir98_AP
  |   |   |   | DWI_dir98_PA
  |   |   |   | SpinEchoFieldMap_PA
  |   |   |   | SpinEchoFieldMap_AP
  |   |   |   | rfMRI_REST_AP
  |   |   |   | rfMRI_REST_PA
  |   |   |   | DWI_dir99_AP
  |   |   |   | DWI_dir99_PA
```