

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\Localizer
 TA:9.2 s PAT:Off Voxel size:1.2x1.2x5.0 mm Rel. SNR:1.00 :f1

Properties

Prio Recon On
 Before measurement
 After measurement
 Load to viewer Off
 Inline movie Off
 Auto store images On
 Load to stamp segments On
 Load images to graphic segments On
 Auto open inline display Off
 Wait for user to start On
 Start measurements single

Routine

Nr. of slice groups 3
 Slices 1
 Dist. factor 20 %
 Position L0.0 A45.0 H0.0 mm
 Orientation Transversal
 Phase enc. dir. A >> P
 AutoAlign ---
 Phase oversampling 0 %
 FoV read 300 mm
 FoV phase 100.0 %
 Slice thickness 5.0 mm
 TR 40.0 ms
 TE 3.00 ms
 Averages 1
 Concatenations 1
 Filter Prescan Normalize, Elliptical filter
 Coil elements HEA;HEP

Contrast

MTC Off
 Magn. preparation None
 Flip angle 15 deg
 Fat suppr. None
 Water suppr. None
 SWI Off
 Averaging mode Short term
 Measurements 1
 Reconstruction Magnitude
 Multiple series Off

Resolution

Base resolution 256
 Phase resolution 75 %
 Phase partial Fourier Off
 Interpolation Off
 PAT mode None
 Image Filter Off
 Distortion Corr. Off
 Unfiltered images Off
 Prescan Normalize On
 Normalize Off
 B1 filter Off

Raw filter	Off
Elliptical filter	On
Mode	Inplane
Geometry	
Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	19.740 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Segments	1
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline		
	Distortion correction	Off
Sequence		
	Introduction	On
	Dimension	2D
	Phase stabilisation	On
	Averaging mode	Short term
	Multi-slice mode	Interleaved
	Asymmetric echo	Allowed
	Contrasts	1
	Bandwidth	260 Hz/Px
	Flow comp.	No
	Allowed delay	0 s
	RF pulse type	Normal
	Gradient mode	Fast
	Excitation	Slice-sel.
	RF spoiling	On
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HEA;HEP
	Acquisition duration	0 ms
	Mode	Off
BOLD		
	Subtract	Off
	Liver registration	Off
	Save images	On
	Autoscaling	Off
	Scaling factor	1
	Offset	0
	Subtrahend	1
	Subtraction indices	
	StdDev	Off
	Std-Dev-Sag	Off
	Std-Dev-Cor	Off
	Std-Dev-Tra	Off
	Std-Dev-Time	Off
	MIP-Sag	Off
	MIP-Cor	Off
	MIP-Tra	Off
	MIP-Time	Off
	Radial MIP	Off
	Save original images	On
	Distortion Corr.	Off
	Contrasts	1
	Save original images	On
	Wash - In	Off
	Wash - Out	Off
	TTP	Off
	PEI	Off
	MIP - time	Off
	Number of radial views	1
	Axis of radial views	L-R
	MPR Sag	Off
	MPR Cor	Off
	MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 Properties

Prio Recon	On
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP
AutoAlign	Head

Contrast

Flip angle	8 deg
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude

Resolution

Base resolution	160
Phase resolution	100 %
Phase partial Fourier	6/8
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	69 %
Slice partial Fourier	6/8

Geometry

Nr. of slab groups	1
--------------------	---

Slabs	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
Multi-slice mode	Sequential
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	24.593 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Weak
Contrasts	1
Bandwidth	540 Hz/Px
RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.

RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Time to center	6.2 s
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\Localizer_aligned
 TA:0:22 PAT:Off Voxel size:1.2x1.2x5.0 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	On
Before measurement	
After measurement	
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal

Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	104.0 ms
TE	3.00 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	15 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane
Geometry	
Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	19.740 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Segments	1
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	2D
Phase stabilisation	On
Averaging mode	Short term
Multi-slice mode	Interleaved
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off

Liver registration	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\BIAS_BC
 TA:0:26 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	250.0 ms
TE	1.01 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BC
Contrast	
Magn. preparation	None
Flip angle	3 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement
Resolution	
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	On
HEP	Off
HEA	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	27.667 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	BC
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1

Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\BIAS_32CH
 TA:0:26 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	250.0 ms
TE	1.01 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HEA;HEP

Contrast

Magn. preparation	None
Flip angle	3 deg

Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement
Resolution	
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	27.667 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On

Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\SpinEchoFieldMap_RL
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None

Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\SpinEchoFieldMap_LR
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms

Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\rfMRI_REST_RL
TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off

Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00

Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\rfMRI_REST_LR
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off

Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px

Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\T1w_MPR
 TA:6:38 PAT:2 Voxel size:0.8x0.8x0.8 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	23.1 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	2400.0 ms
TE	2.12 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize

Coil elements	HEA;HEP
Contrast	
Magn. preparation	Non-sel. IR
TI	1000 ms
Flip angle	8 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement
Resolution	
Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	23.1 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
! Position	L0.0 P3.0 H6.0 mm
! Rotation	90.00 deg
! A >> P	208 mm
! R >> L	180 mm
! F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
ExcitWEns 0 1H	36.889 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	5.000
Physio	
1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1000 ms
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	210 Hz/Px
Flow comp.	No
Echo spacing	7.5 ms
Turbo factor	256
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	

StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\T2w_SPC
 TA:5:57 PAT:2 Voxel size:0.8x0.8x0.8 mm Rel. SNR:1.00 :spc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	3200 ms
TE	563.0 ms
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	None
Water suppr.	None
Restore magn.	Off

Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement
Resolution	
Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Allowed
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off
Geometry	
Nr. of slab groups	1
Slabs	1
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	208
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Table position	P
Restore magn.	Off
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
! Position	L0.0 P3.0 H6.0 mm
! Rotation	90.00 deg
! A >> P	208 mm
! R >> L	180 mm
! F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	138.333 V
! Gain	High
Table position	0 mm
Img. Scale. Cor.	5.000
Physio	
1st Signal/Mode	None
Trigger delay	0 ms
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Bandwidth	744 Hz/Px
Flow comp.	No
Allowed delay	0 s
Echo spacing	3.43 ms
Adiabatic-mode	Off
Turbo factor	314
Echo train duration	1074
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Organ under exam.	None
BOLD	
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off

MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\SpinEchoFieldMap_RL
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %

Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None

Inline		
	Distortion correction	Off
Sequence		
	Introduction	Off
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Bandwidth	2290 Hz/Px
	Echo spacing	0.58 ms
	EPI factor	90
	RF pulse type	Normal
	Gradient mode	Fast
	Triggering scheme	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HEA;HEP
	Acquisition duration	0 ms
BOLD		
	GLM Statistics	Off
	Dynamic t-maps	Off
	Starting ignore meas	0
	Ignore after transition	0
	Model transition states	Off
	Temp. highpass filter	Off
	Threshold	4.00
	Paradigm size	3
	Motion correction	Off
	Spatial filter	Off
	Delay in TR	0 ms
	Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\SpinEchoFieldMap_LR
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties		
	Prio Recon	Off
	Before measurement	
	After measurement	
	Load to viewer	On
	Inline movie	Off
	Auto store images	On
	Load to stamp segments	Off
	Load images to graphic segments	Off
	Auto open inline display	Off
	Wait for user to start	Off
	Start measurements	single
Routine		
	Nr. of slice groups	1
	Slices	72
	Dist. factor	0 %
	Position	L0.0 P3.0 H6.0 mm
	Orientation	T > C-20.0
	Phase enc. dir.	R >> L
	AutoAlign	Head > Brain
	Phase oversampling	0 %
	FoV read	208 mm
	FoV phase	86.5 %

Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\rfMRI_REST_RL
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	

After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None

Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\rfMRI_REST_LR
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None

Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off

Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\tfMRI_WM_RL
 TA:5:01 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8

Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	405
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V

Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\tfMRI_WM_LR
 TA:5:01 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off

Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	405
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On

Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off

Delay in TR 0 ms
Distortion Corr. Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\tfMRI_EMOTION_RL
TA:2:33 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon Off
Before measurement
After measurement
Load to viewer On
Inline movie Off
Auto store images On
Load to stamp segments Off
Load images to graphic segments Off
Auto open inline display Off
Wait for user to start On
Start measurements single

Routine

Nr. of slice groups 1
Slices 72
Dist. factor 0 %
Position L0.0 P3.0 H6.0 mm
Orientation T > C-20.0
Phase enc. dir. R >> L
AutoAlign Head > Brain
Phase oversampling 0 %
FoV read 208 mm
FoV phase 86.5 %
Slice thickness 2.00 mm
TR 720 ms
TE 33.20 ms
Multi-band accel. factor 8
Filter None
Coil elements HEA;HEP

Contrast

MTC Off
Magn. preparation None
Flip angle 52 deg
Fat suppr. Fat sat.
Averaging mode Long term
Measurements 199
Delay in TR 0 ms
Reconstruction Magnitude
Multiple series Off

Resolution

Base resolution 104
Phase resolution 100 %
Phase partial Fourier Off
Interpolation Off
PAT mode None
Distortion Corr. Off
Hamming Off
Prescan Normalize Off
Raw filter Off
Elliptical filter Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast

Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\tfMRI_EMOTION_LR
 TA:2:33 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg

Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	199
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz

Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionA\T1w_vNav_3e
 TA:9:40 PAT:2 Voxel size:0.8x0.8x0.8 mm Rel. SNR:1.00 :tfl

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slab groups	1
Slabs	1

Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	0.80 mm
TR	2400.0 ms
TE 1	1.66 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP
Contrast	
Magn. preparation	Non-sel. IR
TI	1300 ms
Flip angle 1	8.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement
Resolution	
Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Base size phase	256 mm

Base size read	256 mm
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
! Position	L0.0 P3.0 H6.0 mm
! Rotation	90.00 deg
! A >> P	208 mm
! R >> L	180 mm
! F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
ExcitWEns 0 1H	36.889 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	5.000
Physio	
Magn. preparation	Non-sel. IR
TI	1300 ms
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Off
Contrasts	3
Bandwidth 1	740 Hz/Px
Flow comp. 1	No
Echo spacing	8.8 ms
Turbo factor	208
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.

RF spoiling	On
Readout polarity	Positive
Apply moco to	parent and nav
Remeasure	60 TRs
Feedback Delay	108 ms
Moco Ref. Image	Use Temp Ref
Add. grad time	0.0 ms
Apply freq to	parent and nav
Averaging	RMS
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Contrasts	3
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\Localizer
 TA:9.2 s PAT:Off Voxel size:1.2x1.2x5.0 mm Rel. SNR:1.00 :f1

Properties

Prio Recon	On
Before measurement	
After measurement	
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	On

Start measurements	single
Routine	
Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	40.0 ms
TE	3.00 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	15 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane
Geometry	
Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None

Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	19.740 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Segments	1
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	2D
Phase stabilisation	On
Averaging mode	Short term
Multi-slice mode	Interleaved
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz

TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Liver registration	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\AAHScout
 TA:0:14 PAT:3 Voxel size:1.6x1.6x1.6 mm Rel. SNR:1.00 :f1

Properties

Prio Recon	On
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1

Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP
AutoAlign	Head
Contrast	
Flip angle	8 deg
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Resolution	
Base resolution	160
Phase resolution	100 %
Phase partial Fourier	6/8
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	69 %
Slice partial Fourier	6/8
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A45.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
Multi-slice mode	Sequential
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	24.593 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Weak
Contrasts	1
Bandwidth	540 Hz/Px
RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Time to center	6.2 s
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off

Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\Localizer_aligned
 TA:0:22 PAT:Off Voxel size:1.2x1.2x5.0 mm Rel. SNR:1.00 :f1

Properties

Prio Recon	On
Before measurement	
After measurement	
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	104.0 ms
TE	3.00 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	15 deg
Fat suppr.	None
Water suppr.	None

SWI	Off
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	256
Phase resolution	75 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane
Geometry	
Nr. of slice groups	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V

Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	19.740 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Segments	1
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	2D
Phase stabilisation	On
Averaging mode	Short term
Multi-slice mode	Interleaved
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Liver registration	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On

Distortion Corr.	Off
Contrasts	1
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\BIAS_BC
 TA:0:26 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	250.0 ms
TE	1.01 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BC

Contrast

Magn. preparation	None
Flip angle	3 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution	
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	On
HEP	Off
HEA	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	27.667 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	BC
Acquisition duration	0 ms
Mode	Off
BOLD	
Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off

MPR Tra

Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\BIAS_32CH
TA:0:26 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon Off
Before measurement
After measurement
Load to viewer On
Inline movie Off
Auto store images On
Load to stamp segments Off
Load images to graphic segments Off
Auto open inline display Off
Wait for user to start Off
Start measurements single

Routine

Nr. of slab groups 1
Slabs 1
Dist. factor 50 %
Position Isocenter
Orientation Sagittal
Phase enc. dir. A >> P
AutoAlign Head > Brain
Phase oversampling 0 %
Slice oversampling 18.2 %
FoV read 256 mm
FoV phase 100.0 %
Slice thickness 2.00 mm
TR 250.0 ms
TE 1.01 ms
Averages 1
Concatenations 1
Filter None
Coil elements HEA;HEP

Contrast

Magn. preparation None
Flip angle 3 deg
Fat suppr. None
Water suppr. None
Averaging mode Long term
Measurements 1
Reconstruction Magnitude
Multiple series Each measurement

Resolution

Base resolution 128
Phase resolution 100 %
Phase partial Fourier 6/8
Interpolation Off
PAT mode None
Image Filter Off
Distortion Corr. Off
Prescan Normalize Off
Normalize Off
B1 filter Off
Raw filter Off

Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	6/8
Geometry	
Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	88
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.253687 MHz
Correction factor	1
SRFExcit 1H	27.667 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Dark blood	Off
Resp. control	Off
Inline	

Distortion correction	Off
Sequence	
Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	540 Hz/Px
Flow comp.	No
Echo spacing	3 ms
Turbo factor	78
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
Save images	On
Autoscaling	Off
Scaling factor	1
Offset	0
Subtrahend	1
Subtraction indices	
StdDev	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Radial MIP	Off
Save original images	On
Distortion Corr.	Off
Save original images	On
Number of radial views	1
Axis of radial views	L-R
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\SpinEchoFieldMap_RL
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	

Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None

Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0

Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\SpinEchoFieldMap_LR
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off

Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	

Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\rfMRI_REST_LR
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm

TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\rfMRI_REST_RL
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off

Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00

Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\DWI_dir79_RL
 TA:5:11 PAT:Off Voxel size:1.5x1.5x1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	85.7 %
Slice thickness	1.50 mm
TR	3730 ms
TE	76.60 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off

Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off
Geometry	
Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	210 mm
R >> L	180 mm
F >> H	140 mm
Frequency 1H	123.253687 MHz
Correction factor	1
ExtExciterRF 1H	85.623 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On

Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Echo spacing	0.69 ms
EPI factor	120
Gradient mode	Fast
Online multi-band recon.	Remote
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	1
b-value	3000 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off
Calculated bValue	1400 s/mm ²

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\DWI_dir79_LR
 TA:5:11 PAT:Off Voxel size:1.5x1.5x1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	85.7 %
Slice thickness	1.50 mm

TR	3730 ms
TE	76.60 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off
Geometry	
Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off

Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	210 mm
R >> L	180 mm
F >> H	140 mm
Frequency 1H	123.253687 MHz
Correction factor	1
ExtExciteRF 1H	85.623 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Echo spacing	0.69 ms
EPI factor	120
Gradient mode	Fast
Online multi-band recon.	Remote
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	1
b-value	3000 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off
Calculated bValue	1400 s/mm ²

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\DWI_dir84_RL
 TA:5:30 PAT:Off Voxel size:1.5x1.5x1.5 mm Rel. SNR:1.00 :epse

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	85.7 %
Slice thickness	1.50 mm
TR	3730 ms
TE	76.60 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term
Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off
Geometry	
Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved

Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	210 mm
R >> L	180 mm
F >> H	140 mm
Frequency 1H	123.253687 MHz
Correction factor	1
ExtExciteRF 1H	85.623 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Echo spacing	0.69 ms
EPI factor	120
Gradient mode	Fast
Online multi-band recon.	Remote
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	1
b-value	3000 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off
Calculated bValue	1400 s/mm ²

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\DWI_dir84_LR
 TA:5:30 PAT:Off Voxel size:1.5x1.5x1.5 mm Rel. SNR:1.00 :epse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	210 mm
FoV phase	85.7 %
Slice thickness	1.50 mm
TR	3730 ms
TE	76.60 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	78 deg
Fat suppr.	None
Grad. rev. fat suppr.	Enabled
Averaging mode	Long term

Measurements	1
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	140
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Dynamic Field Corr.	Off
Geometry	
Nr. of slice groups	1
Slices	93
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Grad. rev. fat suppr.	Enabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	210 mm
R >> L	180 mm
F >> H	140 mm
Frequency 1H	123.253687 MHz
Correction factor	1

ExtExciterRF LH	85.623 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Resp. control	Off
Inline	
Distortion correction	Off
Sequence	
Introduction	On
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	1700 Hz/Px
Echo spacing	0.69 ms
EPI factor	120
Gradient mode	Fast
Online multi-band recon.	Remote
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
Delay in TR	0 ms
Diffusion mode	Free
Diff. weightings	1
b-value	3000 s/mm ²
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	On
Tensor	Off
Distortion Corr.	Off
b-Value >=	0 s/mm ²
Exponential ADC Maps	Off
Invert Gray Scale	Off
Calculated Image	Off
Calculated bValue	1400 s/mm ²

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\SpinEchoFieldMap_RL
 TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	7080 ms
TE	58.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Grad. rev. fat suppr.	Disabled
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms

Distortion Corr.

Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

\\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\SpinEchoFieldMap_LR
TA:0:28 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epse

Properties

Prio Recon Off
Before measurement
After measurement
Load to viewer On
Inline movie Off
Auto store images On
Load to stamp segments Off
Load images to graphic segments Off
Auto open inline display Off
Wait for user to start Off
Start measurements single

Routine

Nr. of slice groups 1
Slices 72
Dist. factor 0 %
Position L0.0 P3.0 H6.0 mm
Orientation T > C-20.0
Phase enc. dir. R >> L
AutoAlign Head > Brain
Phase oversampling 0 %
FoV read 208 mm
FoV phase 86.5 %
Slice thickness 2.00 mm
TR 7080 ms
TE 58.0 ms
Multi-band accel. factor 1
Filter None
Coil elements HEA;HEP

Contrast

MTC Off
Magn. preparation None
Flip angle 90 deg
Fat suppr. Fat sat.
Grad. rev. fat suppr. Disabled
Averaging mode Long term
Measurements 3
Delay in TR 0 ms
Reconstruction Magnitude
Multiple series Off

Resolution

Base resolution 104
Phase resolution 100 %
Phase partial Fourier Off
Interpolation Off
PAT mode None
Distortion Corr. Off
Hamming Off
Prescan Normalize Off
Raw filter Off
Elliptical filter Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Grad. rev. fat suppr.	Disabled
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
AddCSaCSatNS 1H	40.921 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90

RF pulse type	Normal
Gradient mode	Fast
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\rfMRI_REST_LR
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None

Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm

Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\rfMRI_REST_RL
 TA:5:12 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1

Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	420
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	274
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L

Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\tfMRI_SOCIAL_LR
 TA:3:27 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	274
Delay in TR	0 ms
Reconstruction	Magnitude

Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio		
1st Signal/Mode		None
Magn. preparation		None
Inline		
Distortion correction		Off
Sequence		
Introduction		Off
Averaging mode		Long term
Multi-slice mode		Interleaved
Bandwidth		2290 Hz/Px
Echo spacing		0.58 ms
EPI factor		90
Gradient mode		Fast
Online multi-band recon.		Remote
Triggering scheme		Standard
TX/RX delta frequency		0 Hz
TX Nucleus		None
TX delta frequency		0 Hz
Coil elements		HEA;HEP
Acquisition duration		0 ms
BOLD		
GLM Statistics		Off
Dynamic t-maps		Off
Starting ignore meas		0
Ignore after transition		0
Model transition states		On
Temp. highpass filter		On
Threshold		4.00
Paradigm size		3
Motion correction		Off
Spatial filter		Off
Delay in TR		0 ms
Distortion Corr.		Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\tfMRI_GAMBLING_RL
 TA:3:11 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties		
Prio Recon		Off
Before measurement		
After measurement		
Load to viewer		On
Inline movie		Off
Auto store images		On
Load to stamp segments		Off
Load images to graphic segments		Off
Auto open inline display		Off
Wait for user to start		On
Start measurements		single
Routine		
Nr. of slice groups		1
Slices		72
Dist. factor		0 %
Position		L0.0 P3.0 H6.0 mm
Orientation		T > C-20.0
Phase enc. dir.		R >> L
AutoAlign		Head > Brain

Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	253
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

 \\USER\Lifespan\LS_Phase1b_14-55yo\SessionB\tfMRI_GAMBLING_LR
 TA:3:11 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties	
Prio Recon	Off

Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single
Routine	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Orientation	T > C-20.0
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	208 mm
FoV phase	86.5 %
Slice thickness	2.00 mm
TR	720 ms
TE	33.20 ms
Multi-band accel. factor	8
Filter	None
Coil elements	HEA;HEP
Contrast	
MTC	Off
Magn. preparation	None
Flip angle	52 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	253
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off
Resolution	
Base resolution	104
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Geometry	
Nr. of slice groups	1
Slices	72
Dist. factor	0 %
Position	L0.0 P3.0 H6.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.

Special sat.	None
Special sat.	None
Table position	P
System	
Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.0 P3.0 H6.0 mm
Rotation	90.00 deg
A >> P	208 mm
R >> L	180 mm
F >> H	144 mm
Frequency 1H	123.253687 MHz
Correction factor	1
MBExc 1H	280.027 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000
Physio	
1st Signal/Mode	None
Magn. preparation	None
Inline	
Distortion correction	Off
Sequence	
Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Bandwidth	2290 Hz/Px
Echo spacing	0.58 ms
EPI factor	90
Gradient mode	Fast
Online multi-band recon.	Remote
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
BOLD	
GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0

Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off

SIEMENS MAGNETOM ConnectomS syngo MR D11

Table Of Contents

\\USER

Lifespan	LS_Phase1b_14-55yo	SessionA	Localizer
			AAHScout
			Localizer_aligned
			BIAS_BC
			BIAS_32CH
			SpinEchoFieldMap_RL
			SpinEchoFieldMap_LR
			rfMRI_REST_RL
			rfMRI_REST_LR
			T1w_MPR
			T2w_SPC
			SpinEchoFieldMap_RL
			SpinEchoFieldMap_LR
			rfMRI_REST_RL
			rfMRI_REST_LR
			tfMRI_WM_RL
			tfMRI_WM_LR
			tfMRI_EMOTION_RL
			tfMRI_EMOTION_LR
			T1w_vNav_3e
		SessionB	Localizer
			AAHScout
			Localizer_aligned
			BIAS_BC
			BIAS_32CH
			SpinEchoFieldMap_RL
			SpinEchoFieldMap_LR
			rfMRI_REST_LR
			rfMRI_REST_RL
			DWI_dir79_RL
			DWI_dir79_LR
			DWI_dir84_RL
			DWI_dir84_LR
			SpinEchoFieldMap_RL
			SpinEchoFieldMap_LR
			rfMRI_REST_LR
			rfMRI_REST_RL
			tfMRI_SOCIAL_RL

			tfMRI_SOCIAL_LR
			tfMRI_GAMBLING_RL
			tfMRI_GAMBLING_LR