

Object Edit Tools Insert Help



1	---	Notes ---	
2		Localizer	00:09
3		AAHScout	00:14 AutoAlign Scout
4		Localizer_aligned	00:21
5		BIAS_BC	00:26
6		BIAS_32CH	00:26
7		SpinEchoFieldMap_AP	00:32
8		SpinEchoFieldMap_PA	00:32
9		rfMRI_REST_AP	05:46
10		rfMRI_REST_PA	05:46
11		T1w_MPR	06:38
12		T2w_SPC	05:57
13		dMRI_dir98_AP	05:37
14		dMRI_dir98_PA	05:37
15		SpinEchoFieldMap_AP	00:32
16		SpinEchoFieldMap_PA	00:32
17		rfMRI_REST_AP	05:46
18		rfMRI_REST_PA	05:46
19		dMRI_dir99_AP	05:41
20		dMRI_dir99_PA	05:41



1	---	Notes ---	
2		Local	---
3		AAHS	1) The following is a protocol for a Prisma that is consistent with the Lifespan piloting conducted as part of the HCP at UMN/CMRR.
4		Local	2) The "BIAS*" scans are not currently used in the HCP Pipeline processing, and are therefore "optional" in that regard, although uses may potentially be developed for them in the future.
5		BIAS	3) Note that phase encoding polarity swapping is accomplished via the "Invert RO/PE polarity" flag on the Sequence:Special tab.
6		BIAS	4) Note that the T1/T2w structurals are set up to use an "Adjust Volume" that matches the size and positioning of the SpinEchoFieldMap and rfMRI_REST scans, but for that to work you need to select "Manual" on the Scan Assistant pop-up that will occur when you open these scans.
7		SpinE	3/7/2016
8		SpinE	
9		rfMRI	
10		rfMRI_REST_PA	05:46
11		T1w_MPR	06:38
12		T2w_SPC	05:57
13		dMRI_dir98_AP	05:37
14		dMRI_dir98_PA	05:37
15		SpinEchoFieldMap_AP	00:32
16		SpinEchoFieldMap_PA	00:32
17		rfMRI_REST_AP	05:46
18		rfMRI_REST_PA	05:46
19		dMRI_dir99_AP	05:41
20		dMRI_dir99_PA	05:41

Object Edit Tools Insert Help



1	---	Notes ---	
2		Localizer	00:09
3		AAHScout	00:14 AutoAlign Scout
4		Localizer_aligned	00:21
5		BIAS_BC	00:26
6		BIAS_32CH	00:26
7		SpinEchoFieldMap_AP	00:32
8		SpinEchoFieldMap_AP	
9		rfMRI	
10		rfMRI_REST_PA	05:46
11		T1w_MPR	06:38
12		T2w_SPC	05:57
13		dMRI_dir98_AP	05:37
14		dMRI_dir98_PA	05:37
15		SpinEchoFieldMap_AP	00:32
16		SpinEchoFieldMap_PA	00:32
17		rfMRI_REST_AP	05:46
18		rfMRI_REST_PA	05:46
19		dMRI_dir99_AP	05:41
20		dMRI_dir99_PA	05:41

SpinEchoFieldMap_AP
Use "Free echo spacing" option on Sequence tab to match
echo spacing to BOLD scans.

Object Edit Tools Insert Help



1	---	Notes ---	
2		Localizer	00:09
3		AAHScout	00:14 AutoAlign Scout
4		Localizer_aligned	00:21
5		BIAS_BC	00:26
6		BIAS_32CH	00:26
7		SpinEchoFieldMap_AP	00:32
8		SpinEchoFieldMap_PA	00:32
9		rfMRI_REST_AP	05:46
10		rfMRI_REST_PA	
11		No PF (partial Fourier) and TR=800 ms (rather than 720 ms), which are small differences relative to the Lifespan piloting on the UMN/CMRR Prisma.	
12		T2w_SPC	05:57
13		dMRI_dir98_AP	05:37
14		dMRI_dir98_PA	05:37
15		SpinEchoFieldMap_AP	00:32
16		SpinEchoFieldMap_PA	00:32
17		rfMRI_REST_AP	05:46
18		rfMRI_REST_PA	05:46
19		dMRI_dir99_AP	05:41
20		dMRI_dir99_PA	05:41