


## Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Slice group 1   FoV read 210   mmSlices 92  FoV phase 100.0   %Dist. factor 0   %Slice thickness 1.50   mmPosition L0.0 P3.0 H6.0  TR 3230   msOrientation T > C-20.0  TE 89.20   msPhase enc. dir. A >> P  AutoAlign Head > Brain  Phase oversampling 0   %Multi-band accel. factor 4  

Filter None

Coil elements HEA,HEP



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK


Cancel

Virtual Coils...

Help

## Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Slice group 1   FoV read 210   mmSlices 92  FoV phase 100.0   %Dist. factor 0   %Slice thickness 1.50   mmPosition L0.0 P3.0 H6.0  TR 3230   msOrientation T > C-20.0  TE 89.20   msPhase enc. dir. A >> P  AutoAlign Head > Brain  Phase oversampling 0   %Multi-band accel. factor 4  Filter None Coil elements HEA,HEP 

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

## Common

## Dynamic

TR 3230 ms

Fat suppr. None


TE 89.20 ms

Grad. rev. fat suppr. Enabled

MTC ☐

Magn. preparation None

Flip angle 78 deg



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

## Common

## Dynamic

TR 3230 ms

Fat suppr. None

TE 89.20 ms

Grad. rev. fat suppr. Enabled

MTC ☐

Magn. preparation None

Refocus flip angle 180 deg

Use the toggle selector to get to  
"Refocus flip angle".

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

**Common****Dynamic**

Averaging mode Long term ▼

Measurements 1 ▼

Delay in TR 0 ▼ ms

Reconstruction Magnitude ▼

Multiple series Off ▼

Routine

**Contrast**

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Common

iPAT

Filter Image

Filter Rawdata

FoV read 210 mm


FoV phase 100.0 %

Slice thickness 1.50 mm

Base resolution 140

Phase resolution 100 %

Phase partial Fourier 6/8

Interpolation ☐

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

\\USER\\head\\Harms\\CCF\_PrismalMRI\_dir98\_AP

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Common

**iPAT**

Filter Image

Filter Rawdata

PAT mode

None



Routine

Contrast

**Resolution**

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Common

iPAT

**Filter Image**

Filter Rawdata

Prescan Normalize



Distortion Corr.



Dynamic Field Corr.



Routine

Contrast

**Resolution**

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help



**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Common

iPAT

Filter Image

Filter Rawdata

Raw filter



Elliptical filter



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

## Common

## Saturation

## Navigator

Slice group 1

FoV read 210 mm

Slices 92

FoV phase 100.0 %

Dist. factor 0 %

Slice thickness 1.50 mm

Position L0.0 P3.0 H6.0

TR 3230 ms

Orientation T &gt; C-20.0


Phase enc. dir. A &gt;&gt; P

Multi-slice mode Interleaved

Series Interleaved

Phase oversampling 0 %

Multi-band accel. factor 4



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

**Common****Saturation****Navigator**

Fat suppr. None



Grad. rev. fat suppr. Enabled



Sat. region



Special sat. None



Routine

Contrast

Resolution

**Geometry**

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

\\USER\\head\\Harms\\CCF\_Prismal\\dMRI\_dir98\_AP

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Common

Saturation

**Navigator**

Navigator



Routine

Contrast

Resolution

**Geometry**

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Coils

Miscellaneous

Adjustments

Adjust Volume

pTx Volumes

Tx/Rx



HEA

HEP

Body

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Coils

**Miscellaneous**

Adjustments

Adjust Volume

pTx Volumes

Tx/Rx

Coil Combine Mode **Sum of Squares**Positioning mode **REF**Table position **H** **0** mmMatrix Optimization **Off**Coil Focus **Flat****Image Numbering**MSMA **S - C - T**Sagittal **R >> L**Coronal **A >> P**Transversal **F >> H**AutoAlign **Head > Brain**Coil Select Mode **Off - All**

Routine

Contrast

Resolution

Geometry

**System**

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

**Coils****Miscellaneous****Adjustments****Adjust Volume****pTx Volumes****Tx/Rx**B0 Shim mode **Standard** ▼B1 Shim mode **TrueForm** ▼Adjustment Tolerance **Auto** ▼Adjust with body coil ☐Confirm freq. adjustment ☐Assume Dominant Fat ☐Assume Silicone ☐

Tx Ref [Nucleus]    Ref.

? Ref. amplitude 1H	0.000

**Reset****Routine****Contrast****Resolution****Geometry****System****Physio****Diff****Sequence****OK****Cancel****Virtual Coils...****Help**

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

**Coils****Miscellaneous****Adjustments****Adjust Volume**

pTx Volumes

Tx/Rx

Position L0.0 P3.0 H6.0



Orientation T &gt; C-20.0



Rotation 0.00 deg

R &gt;&gt; L 210 mm

A &gt;&gt; P 210 mm

F &gt;&gt; H 138 mm

Reset

**Routine****Contrast****Resolution****Geometry****System****Physio****Diff****Sequence**

OK

Cancel

Virtual Coils...

Help



\\USER\\head\\Harms\\CCF\_PrismalMRI\_dir98\_AP

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Coils

Miscellaneous

Adjustments

Adjust Volume

**pTx Volumes**

Tx/Rx

B1 Shim mode

TrueForm



pTx Volume



+

-

Routine

Contrast

Resolution

Geometry

**System**

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37 PM: REF PAT: Off Voxel size: 1.5×1.5×1.5 mm Rel. SNR: 1.00 : epse

Coils

Miscellaneous

Adjustments

Adjust Volume

pTx Volumes

Tx/Rx

## Transmitter

## Receiver

Frequency 1H 123.257575 MHz

Gain High

? Ref. amplitude 1H 0.000 V

Correction factor 1

Img. Scale. Cor. 1.000

Puls Amplitude V


ExtExciteRF 1H	114.337
ExtRefocRF 1H	126.069
MBExc 1H	457.346
MBRef 1H	504.278

Make sure that "MBExc 1H" and "MBRef 1H" are not maxed out, otherwise you will get clipping of your RF pulses (although, as of the R013 release, a warning message should be generated at run time if this is the case).

Power issues are controlled by the values/settings of "Excite/refocus pulse duration", "MB RF phase scramble", and "Time-shifted MB RF" on the Sequence:Special tab.

Reset

MBExc 1H



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37 PM: REF PAT: Off Voxel size: 1.5×1.5×1.5 mm Rel. SNR: 1.00 : epse

Coils

Miscellaneous

Adjustments

Adjust Volume

pTx Volumes

Tx/Rx

## Transmitter

## Receiver

Frequency 1H 123.257575 MHz

Gain High

? Ref. amplitude 1H 0.000 V

Correction factor 1


Img. Scale. Cor. 1.000

Puls Amplitude V

ExtExciteRF 1H	114.337
ExtRefocRF 1H	126.069
MBExc 1H	457.346
MBRef 1H	504.278

Reset

MBRef 1H



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

**Signal1**

PACE

1st Signal/Mode

None

TR 3230 ms

Multi-band accel. factor

4

Routine

Contrast

Resolution

Geometry

System

**Physio**

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Signal1

**PACE**


Resp. control

Off

Multi-band accel. factor

4

Multi-band accel. factor



Routine

Contrast

Resolution

Geometry

System

**Physio**

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Neuro

Body

Diffusion mode **Free** ...

Diff. directions 98 ...

Diffusion Scheme Monopolar

Diff. weightings 2

Diff. weighted images ☒Trace weighted images ☐ADC maps ☐FA maps ☐Mosaic ☒Tensor ☐

b-value

b-value 1	0
b-value 2	3000

Note: "b-value 2" = 3000 reflects the \*maximum\* b-value. By scaling the vectors in the gradient table appropriately, other b-values ("shells") are acquired when the scan runs. E.g., for the particular 98 direction file used here, the norm of the gradient direction vectors was either 1.0 or  $\sqrt{1/2}$ , resulting in shells of  $3000 * [1.0 \ 1/2] = \{3000, 1500\}$  (not counting the interspersed b=0 scans).

N.B. Need to set "Diff. weightings"  $\geq 2$  to activate option to select "Mosaic" for the output format. Here, we've set "b-value 1" = 0, which results in a single additional b=0 volume acquired prior to cycling through the specified gradient table (which itself includes interspersed b=0 scans).

Noise level 40

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

Neuro

Body

Diffusion mode **Free**

Diff. directions 98

Diffusion Scheme Monopolar

Diff. weightings 2

b-value

b-value 1	0
b-value 2	3000

Diff. weighted images ☒Trace weighted images ☐ADC maps ☐Exponential ADC Maps ☐FA maps ☐Invert Gray Scale ☐Calculated Image ☐b-Value >= 0 s/mm<sup>2</sup>

Noise level 40

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37    PM: REF    PAT: Off    Voxel size: 1.5×1.5×1.5 mm    Rel. SNR: 1.00    : epse

**Part 1**

Part 2

**Special**

Introduction



Bandwidth 1700 Hz/Px

Averaging mode Long term

Multi-slice mode Interleaved

Free echo spacing

Echo spacing 0.69 ms

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

**Sequence**

OK

Cancel

Virtual Coils...

Help



**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

Special

Gradient mode

Performance

RF spoiling



EPI factor

140

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

Special

Excite pulse duration

3840

us

Single-band images



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

Special

Refocus pulse duration

7680

us



Use toggle to get to "Refocus pulse duration".

Single-band images

Another toggle here reveals other parameters.  
In particular, for dMRI want to use "SENSE1 coil combine".

Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

**Special**

Refocus pulse duration

7680

us



MB LeakBlock kernel



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

**Sequence**

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

**Special**

Refocus pulse duration

7680

us



MB RF phase scramble



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

**Sequence**

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

**Special**

Refocus pulse duration

7680

us



Time-shifted MB RF



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

**Sequence**

OK

Cancel

Virtual Coils...

Help

## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

Special

Refocus pulse duration

7680

us



SENSE1 coil combine



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help

**Details**

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

**Special**

Refocus pulse duration

7680

us



Log physiology to file



Online multi-band recon.

Online



FFT scale factor

1.00



Routine

Contrast

Resolution

Geometry

System

Physio

Diff

**Sequence**

OK

Cancel

Virtual Coils...

Help



## Details

TA: 5:37

PM: REF

PAT: Off

Voxel size: 1.5×1.5×1.5 mm

Rel. SNR: 1.00

: epse

Part 1

Part 2

Special

Refocus pulse duration

7680

us



Invert RO/PE polarity



Online multi-band recon.

Online



FFT scale factor

1.00



N.B. Checking "Invert RO/PE polarity" is the mechanism we use to obtain a matched scan with opposite phase encoding polarity (while leaving "Phase enc. dir." on the Routine tab set to "A >> P").

Routine

Contrast

Resolution

Geometry

System

Physio

Diff

Sequence

OK

Cancel

Virtual Coils...

Help